

Mepilex® Application Guide

Heel

Examples: • Dressing a wound on a heel • Protection for blisters
• Cushioning to prevent pressure ulcers



Between toes

Examples: • Dressing of a wound between toes



Mepilex 10 x 10 cm



Tip of a toe

Examples: • Cushioning out a tip of a toe • Toe nail removal



Mepilex 10 x 10 cm



Fingers and toes

Examples: • Finger injuries • Toe injuries




Mepilex 10 x 10 cm



Hands

Examples: • Hand contraction • Knuckle or palm wounds



Mepilex 10 x 10 or 10 x 20cm



Elbows

Examples: • Cushioning to prevent pressure ulcers
• Post treatment of bursitis



Mepilex 10 x 10




Axilla

Examples: • Lymph nodes • Exuding wounds




Mepilex 20 x 20



Ears

Examples: • Cushioning to prevent pressure ulcers on the ears
• Wounds around the ear



Mepilex 10 x 10




Around tubes and catheters

Examples: • At PEG, suprapubic catheter and drainage sites
• Can also be used around external fixation pins



Mepilex 10 x 10



MEPILEX ASSORTMENT

Art. no	Size cm	Pcs/shelf	Pcs/transp. cont.
294100	10 x 10	5	70
294200	10 x 20	5	45
294300	15 x 15	5	25
294400	20 x 20	2	20
294500	20 x 50	2	12
Mepilex Heel			
288100	13 x 20	5	25